

Amplicon-based targeted sequencing of single circulating tumor cells

RARECYTE

Nolan G. Ericson¹, Vidushi Kapoor², Guoying Liu², Arturo B. Ramirez¹, Alisa C. Clein³, Celestia S. Higano³, Daniel E. Sabath³, Eric P. Kaldjian¹, Tad George¹
¹RareCyte, Inc., Seattle, WA, ²Paragon Genomics, Hayward, CA, ³University of Washington, Seattle, WA



Background

RareCyte has developed technology for visual identification and single cell retrieval of rare cells in blood, including circulating tumor cells (CTCs). There is increasing interest in mutational analysis of circulating tumor cells (CTCs) as a liquid biopsy application. Because of the minuscule amount of DNA present in a single cell (~6 pg), whole genome amplification (WGA) is typically performed prior to next generation sequencing (NGS) library preparation. Existing WGA methods have inherent amplification biases leading to non-uniform genome coverage that can cause dropout of desired targets, as well as elevated error rates that can lead to false positive mutations. Paragon Genomics has developed the CleanPlex[®] amplicon-based target enrichment technology for ultra-high multiplexing of amplicons with high coverage uniformity. Here we apply CleanPlex without pre-amplification for single CTC targeted DNA sequencing.

Methods

Cultured A549 cells as “model CTCs” were spiked into whole blood from healthy donors. Blood was processed by AccuCyte[®] Sample Preparation System onto slides. Multi-parameter immunofluorescence and automated imaging with the CyteFinder[®] Instrument were used to identify CTCs. CTCs were retrieved by CytePicker[®] into PCR tubes and either amplified by WGA (PicoPLEX[®]) or lysed in a PCR-compatible lysis buffer. In clinical samples, patient-matched plasma DNA from RareCyte Blood Collection Tubes or Streck Cell-Free DNA BCT[®] was also isolated. WGA products, cell lysates and plasma DNA were used as template for the CleanPlex OncoZoom[®] Cancer Hotspot Panel; additional PCR cycles were added during target amplification to compensate for low DNA input in the non-WGA single cell samples. Libraries were sequenced on the Illumina MiSeq.

Results

Increased uniformity of coverage with decreased target dropout was observed with CleanPlex NGS library prep direct from cell lysate (non-WGA). Median read depth increased 48-fold when compared to the WGA method. On average, 23 out of 29 (79%) variants present in bulk A549 genomic DNA were observed with the non-WGA method in single mCTCs, while 15 of 29 (52%) were observed after WGA. Additionally, the false positive error frequency of non-WGA samples was 8-fold lower than in the WGA samples. CTCs and cell-free DNA from two metastatic breast cancer patients were also sequenced using OncoZoom. In one patient, PIK3CA E542K, a well-documented oncogenic mutation, was observed in 3 of 5 CTCs. In the other, ERBB2 L755S, a known HER2-reactivating mutation associated with chemoresistance, was found in 3 of 5 CTCs. Both patients' variants were also seen at similar allelic frequencies in RareCyte and Streck cell-free DNA.

Conclusions

Single-cell CleanPlex amplicon-based sequencing without prior WGA resulted in libraries with more complete and consistent coverage and lower error frequencies, enabling efficient and accurate assessment of somatic mutations in CTCs. This approach shows promise for cell-based liquid biopsy diagnostic applications.

Sample preparation and analysis workflow

