

Xiaoyu An¹, Xiaolong Tu², Xinhe Feng², Jianjian Zhang², Dandan Song², Haijuan Yu², Henry Li¹

¹Crown Bioscience Inc., 16550 West Bernardo Drive, Building 5, Suite 525, San Diego, CA 92127;

²Crown Bioscience Taicang, Science & Technology Innovation Park, No.6 Beijing West Road, Taicang City, Jiangsu Province, P.R. China 215400

INTRODUCTION


Cell therapy for cancers has gained great momentum as a promising new treatment approach. In particular, autologous T cells modified with a chimeric antigen receptor (CAR-T) has recently been regulatory-approved for the treatment of patients with CD19+ acute B cell leukemia^{1,2}. However, effective monitoring and characterization of the CAR-T cells after administration, in order to understand the treatment effect and underlying mechanism, remains an urgent unmet need for CAR-T cell therapy in the clinic. Currently, there are two major methods to monitor the presence of CAR-T cells in patients:

1. Flow cytometry (a single cell based assay recognizing the tag engineered in the CAR molecule), but with extremely poor sensitivity ~1%.
2. Real-Time PCR (amplifying the engineered tag sequences), being more sensitive, but without the ability to identify and characterize individual CAR-T cells^{3,4}.

Thus, a method with high sensitivity, and also the ability to characterize individual CAR-T cells, is in high demand.

METHODS

Fig 1. CAR-T sample preparation, staining, and analysis workflow



We are attempting to establish a novel platform, based on rare cell detection in a pathology format, to detect and characterize rare CAR-T cells in patients treated by CAR-T, using RareCyte technology.

First, 1×10^6 white blood cells (PBMC from C57BL/6 mice), spiked with the serial dilution of 0-250 human breast cancer MCF-7 cells, were smeared on a pathology slide and subjected to staining by antibodies against human cytokeratin 18 (Abcam, ab181597), human EpCAM (Abcam, ab8666), and mouse CD45 (R&D, AF114), followed by analysis by either flow cytometry or the RareCyte system, in order to test the detection sensitivity and quantification dynamic ranges.

Second, purified CD19 CAR-T cells, which were mixed into CD3+ T cells, were analyzed by flow cytometry and the RareCyte system individually, and the detection rate was compared.

Lastly, whole blood obtained from a DLBCL patient treated with CD19 CAR-T cells (100 μ L) was serially diluted with whole blood from healthy donors and then subjected to the same two analyses as detailed above with CAR-T specific staining (CD3-FITC, CD19-CAR-PE, DAPI (nucleus)), as controlled by the whole blood from healthy donors.

RESULTS

Fig 2. Representative staining images of the breast cancer cell line MCF-7 spiked-in mouse WBC




Fig 3. Linear regression and sensitivity analysis







Fig 4. Comparison of detection viability between RareCyte and flow cytometry with purified CD19 CAR-T cells mixed with CD3+ T cells



C

| Cell Type | Field | | | | | | | | | | Average |
|-----------------------------|-------|-------|-------|-------|-------|-------|-----|-------|-------|-------|---------|
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | |
| CD3+ Cells | 433 | 461 | 46 | 117 | 285 | 170 | 122 | 231 | 241 | 346 | 245.20 |
| CAR-T Cells | 251 | 262 | 25 | 72 | 157 | 97 | 61 | 124 | 132 | 197 | 137.80 |
| %CAR-T Cells/ CD3+ Cells | 57.97 | 56.83 | 54.35 | 61.54 | 55.09 | 57.06 | 50 | 53.68 | 54.77 | 56.94 | 56.20 |

Fig 5. CAR-T cell identification with RareCyte and flow cytometry in treated patient blood



C

| Detection Platform | CAR-T Treated Patient Blood /Serial Dilution | | | | | | Healthy Blood (Control) |
|--|--|------|------|-------|-------|--------|-------------------------|
| | Undiluted | 1:5 | 1:25 | 1:125 | 1:625 | 1:3125 | |
| RareCyte (CAR-T+ Cell Number) | 141 | 12 | 3 | 0 | 1 | 5 | 0 |
| Flow Cytometry (%CAR-T+ in CD3+ T Cells) | 10.26 | 0.98 | 0.38 | 0.27 | 0.27 | 0.18 | 0.21 |

SUMMARY

- Preliminary data from the spiked study indicates that our method has high sensitivity down to 2 cells (>90% recovery), and a good linear range of quantification. The sensitivity (<1/1-million) is at least 3-orders of magnitude above that of flow cytometry (~1%)
- Our CAR-T treated clinical data demonstrated the detection and quantification of CAR-T cells in patients, whereas flow cytometry has limitations. In addition, a single CAR-T can be picked up for further characterization, e.g. single cell genomics
- Our method for detection and monitoring offers a great opportunity for the characterization of CAR-T cells in the clinical setting, in order to monitor patient treatment and understand treatment mechanisms, that is otherwise unavailable. This could become a useful diagnostic to provide guidance on the treatment and prognosis of patients

REFERENCES

1. Fesnak *et al.* Engineered T cells: the promise and challenges of cancer immunotherapy. *Nat Rev Cancer* 2016;16(9):566-81.
2. Jackson *et al.* Driving CAR T-cells forward. *Nat Rev Clin Oncol* 2016;13(6):370-83.
3. Jena *et al.* Chimeric antigen receptor (CAR)-specific monoclonal antibody to detect CD19-specific T cells in clinical trials. *PLoS One* 2013;8(3):e57838.
4. Wang *et al.* Establishment of a Quantitative Polymerase Chain Reaction Assay for Monitoring Chimeric Antigen Receptor T Cells in Peripheral Blood. *Transplant Proc* 2018;50(1):104-9.